



Artisan Scientific

QUALITY INSTRUMENTATION ... GUARANTEED

Looking for more information?

Visit us on the web at <http://www.artisan-scientific.com> for more information:

- Price Quotations
- Drivers
- Technical Specifications, Manuals and Documentation

Artisan Scientific is Your Source for Quality New and Certified-Used/Pre-owned Equipment

- Tens of Thousands of In-Stock Items
- Hundreds of Manufacturers Supported
- Fast Shipping and Delivery
- Leasing / Monthly Rentals
- Equipment Demos
- Consignment

Service Center Repairs

Experienced Engineers and Technicians on staff in our State-of-the-art Full-Service In-House Service Center Facility

InstraView™ Remote Inspection

Remotely inspect equipment before purchasing with our Innovative InstraView™ website at <http://www.instraview.com>

We buy used equipment! We also offer credit for Buy-Backs and Trade-Ins

Sell your excess, underutilized, and idle used equipment. Contact one of our Customer Service Representatives today!

Talk to a live person: 888-88-SOURCE (888-887-6872) | Contact us by email: sales@artisan-scientific.com | Visit our website: <http://www.artisan-scientific.com>

OVENS (VACUUM)

Lindberg/Blue M Vacuum Ovens

Lindberg/Blue M vacuum ovens are designed for drying, curing, outgassing, aging, process control and other applications which require elevated temperature in reduced atmospheres or vacuum/purge with non-flammable and inert atmospheres. Select from three sizes, all with digital electronic control, built-in overtemperature protection, and a fully flexible vacuum/purge/release system for a range of uses. Single setpoint digital control, front-mounted. Simultaneous digital display of setpoint and actual temperature with push-button view of overtemperature setpoint, resolution to within 1°C. Keypad data entry and pull-out control module for easy service. Vacuum system has a 1" (25.4 mm) manifold, chamber rear, exterior, connects to pump or in-house vacuum source. Includes inert gas injection valve, fresh air inlet, vacuum exhaust port, controllable vacuum release vent port. Heavy gauge steel construction, dual wall, with fiberglass insulation for minimum heat loss in walls and door. Painted exterior or deluxe stainless steel finish. For use with non-flammable, non-corrosive inert gases only. Vacuum gauge, vacuum, gas and inlet valves located on front. High temperature silicone gaskets are supplied standard. Solvent resistant Buna-N high temperature gaskets may be ordered for specific applications where silicone is not preferred and temperatures are limited to +160°C. Viton-B gaskets are also available for applications involving acids. Pumps, cold traps, hoses, vacuum gauges, clamps and fittings must be ordered separately and require customer installation. Lindberg/Blue M. (GM)



O4388

Specifications:

Specifications	O4388-1/-2/-7	O4388-3/-4/-9	O4388-5/-6/-11
Chamber Volume	18.6 Lt	42.5 Lt	127.5 Lt
Chamber (WxHxD")	9 x 14 x 9	12 x 18 x 12	18 x 24 x 18
Temp. Range (° C)	up to 260° C	up to 260° C	up to 260° C
Temp. Uniformity @ 100° C	±3.5 %	±3.5 %	±3.5 %
Rise Time to 200° C	64 min	120 min.	165 min.
Vacuum Capability	10 microns	10 microns	10 microns
Shipping Weight	150 Lbs	250 Lbs	325 Lbs

Cat. No.	Mfg. No.	Description	Qty	Price
O4388-1	VO914A	18.6 Lt Vacuum Oven w/Painted Exterior. 120V	1	\$2,695.00
O4388-2	VO914C	18.6 Lt Vacuum Oven w/Painted Exterior. 208/240V	1	\$2,695.00
O4388-3	VO1218A	42.5 Lt Vacuum Oven w/Painted Exterior. 120V	1	\$4,780.00
O4388-4	VO1218C	42.5 Lt Vacuum Oven w/Painted Exterior. 208/240V	1	\$4,780.00
O4388-5	VO1824A	127.5 Lt Vacuum Oven w/Painted Exterior. 120V	1	\$10,620.00
O4388-6	VO1824C	127.5 Lt Vacuum Oven w/Painted Exterior. 208/240V	1	\$10,620.00
O4388-7	VO914SA	18.6 Lt Vacuum Oven w/ Stainless Steel Exterior. 120V	1	\$4,120.00
O4388-9	VO1218SA	42.5 Lt Vacuum Oven w/ Stainless Steel Exterior. 120V	1	\$6,360.00
O4388-11	VO1824SA	127.5 Lt Vacuum Oven w/ Stainless Steel Exterior. 120V	1	\$12,725.00

Blue Textured Nitrile Gloves, Lightly Powdered. HIGH FIVE

The HIGH FIVE® nitrile glove is 100% Latex-Free, lightly Powdered and contains no natural rubber proteins. This soft and supple glove is odor-free, and offers an unequaled fit and feel compared to conventional exam gloves. Our proprietary low modulus formulation reduces hand fatigue and offers superior puncture resistance and chemical resistance compared to latex or vinyl gloves. No MBT's, Thiriums or Carbamates. 100% inspected and approved for medical use. The HIGH FIVE® powder-free nitrile is the ideal solution for individuals allergic to latex. Hand packed in convenient dispenser boxes; color-coded by size. Other features include: 100% tested for pinholes and defects, ISO 9002 certification, 1.5 AQL pinhole rating, 5-6 mil thickness, non-sterile and ambidextrous, 9 1/2 inches length with a beaded cuff for easing donning. High Five. (CGLJMT)



N811

Cat. No.	Mfg. No.	Description	Size	Color	Qty/bx	Price/bx	Qty/cs	Price/cs
N811	N811	Nitrile, Powder-Free Textured	Small	Blue	100	\$18.29	20 boxes	\$318.00
N812	N812	Nitrile, Powder-Free Textured	Medium	Blue	100	\$18.29	20 boxes	\$318.00
N813	N813	Nitrile, Powder-Free Textured	Large	Blue	100	\$18.29	20 boxes	\$318.00
N814	N814	Nitrile, Powder-Free Textured	X-Large	Blue	100	\$18.29	20 boxes	\$318.00



Artisan Scientific

QUALITY INSTRUMENTATION ... GUARANTEED

Looking for more information?

Visit us on the web at <http://www.artisan-scientific.com> for more information:

- Price Quotations
- Drivers
- Technical Specifications, Manuals and Documentation

Artisan Scientific is Your Source for Quality New and Certified-Used/Pre-owned Equipment

- Tens of Thousands of In-Stock Items
- Hundreds of Manufacturers Supported
- Fast Shipping and Delivery
- Leasing / Monthly Rentals
- Equipment Demos
- Consignment

Service Center Repairs

Experienced Engineers and Technicians on staff in our State-of-the-art Full-Service In-House Service Center Facility

InstraView™ Remote Inspection

Remotely inspect equipment before purchasing with our Innovative InstraView™ website at <http://www.instraview.com>

We buy used equipment! We also offer credit for Buy-Backs and Trade-Ins

Sell your excess, underutilized, and idle used equipment. Contact one of our Customer Service Representatives today!

Talk to a live person: 888-88-SOURCE (888-887-6872) | Contact us by email: sales@artisan-scientific.com | Visit our website: <http://www.artisan-scientific.com>